

Datasheet

KHDRBS1 monoclonal antibody (M03), clone 1A4

Catalog Number: H00010657-M03

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant KHDRBS1.

Clone Name: 1A4

Immunogen: KHDRBS1 (AAH00717, 176 a.a. ~ 275 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

ILGPQGNTIKRLQEETGAKISVLGKGSMRDKAKEEELR
KGGDPKYAHLNMDLHVFIIEVFGPPCEAYALMAHAMEE
VKKFLVPDMMDDICQEQLFLELSYLN

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, S-ELISA, WB-Ce, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10657

Gene Symbol: KHDRBS1

Gene Alias: FLJ34027, Sam68, p62