

## Datasheet

### KDM5B monoclonal antibody (M02A), clone 1G10

**Catalog Number:** H00010765-M02A

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a partial recombinant KDM5B.

**Clone Name:** 1G10

**Immunogen:** KDM5B (NP\_006609, 231 a.a. ~ 319 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

AKRMRAEAMNIKIEPEETTEARTHNLRRRMGCPTPKC  
ENEKEMKSSIKQEPIERKDYIVENEKEKPKSRSKKATN  
AVDLYVCLLCGSGN

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB-Re

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG2a Kappa

**Storage Buffer:** In ascites fluid

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 10765

**Gene Symbol:** KDM5B

**Gene Alias:** FLJ10538, FLJ12459, FLJ12491, FLJ16281, FLJ23670, KDM5B, PLU-1, PUT1, RBBP2H1A