

Datasheet

HLA2 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00011148-B01P

Regulatory Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human HHLA2 protein.

Immunogen: HHLA2 (NP_009003.1, 1 a.a. ~ 414 a.a) full-length human protein.

Sequence:

MKAQTALSFFLILITSLSGSQGIFPLAFFIYVPMNEQIVIG
RLDEDIILPSSFERGSEVVIHWKYQDSYKVHSYYKGS
HLESQDPRYANRTSLFYNEIQNGNASLFFRRVSLLE
GIYTCYVGTAIQVITNKVVLKGVFLTPVMKYEKRNTNS
FLICSVLSVYPRPIITWKMDNTPISENNMEETGSLDSFS
INSPLNITGSNSSYECTIENSLKQWTGRWTMKDGLH
KMQSEHVSLSCQPVNDYFSPNQDFKVTWSRMKSGTF
SVLAYYLSSSQNTIINESRFSWNKELINQSDFSMNLMD
LNLSDSGEYLCNISSDEYTLTIHTVHVEPSQETASHNK
GLWILVPSAILAAFLIWSVKCCRAQLEARRSRHPADG
AQQERCCVPPGERCPSAPDNGEENVPLSGKV

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 11148

Gene Symbol: HHLA2

Gene Alias: -