

Datasheet

STK38 monoclonal antibody (M03), clone 6F1

Catalog Number: H00011329-M03

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant STK38.

Clone Name: 6F1

Immunogen: STK38 (AAH12085, 365 a.a. ~ 465 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

EHRIGAPGVVEIKSNSFFEGVDWEHIRERPAAISIEIKSI
DDTSNFDEFPEKDILKPTVATSNHPETDYKNKDWVFIN
YTYKRFEGLTARGAIPSYMKAAK

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 11329

Gene Symbol: STK38

Gene Alias: NDR, NDR1