

Datasheet

MCF2L monoclonal antibody (M01), clone 1D11

Catalog Number: H00023263-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant MCF2L.

Clone Name: 1D11

Immunogen: MCF2L (AAH20208, 220 a.a. ~ 319 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

GTELAETELPNDVQSTSSVLCAHTEKKDKAKEDLRLAL
KEGHSVLESLRELQAEGSEPSVNQDQLDNQATVQRLL
AQLNETEAAFDEFWAKHQKLEQCL

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 23263

Gene Symbol: MCF2L

Gene Alias: ARHGEF14, DBS, FLJ12122, KIAA0362, OST