

## Datasheet

### CIDEB monoclonal antibody (M01), clone 2A10

**Catalog Number:** H00027141-M01

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a partial recombinant CIDEB.

**Clone Name:** 2A10

**Immunogen:** CIDEB (NP\_055245, 3 a.a. ~ 110 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

YLSALNPSDLLRSVSNISSEFGRRVWTSAPPPQRPFR  
VCDHKRTIRKGLTAATRQELLAKALETLLLNGVLTLVLE  
EDGTAVDSEDFQLEDDTCLMVLQSGQSWSP

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG1 Kappa

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 27141

**Gene Symbol:** CIDEB

**Gene Alias:** -