

## Datasheet

### RPS6KA6 monoclonal antibody (M02), clone 6F2

**Catalog Number:** H00027330-M02

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a partial recombinant RPS6KA6.

**Clone Name:** 6F2

**Immunogen:** RPS6KA6 (NP\_055311.1, 636 a.a. ~ 745 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

RIGNGKFSLSGGNWDNISDGAKDLLSHMLHMDPHQR  
YTAEQILKHSWITHRDQLPNDQPKRNDVSHVVKGAMV  
ATYSALTHKTFQPVLEPVAASSLAQRRSMKKRTSTGL

**Host:** Mouse

**Reactivity:** Human, Rat

**Applications:** ELISA, IHC-P, WB-Ce, WB-Re  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG2a Kappa

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 27330

**Gene Symbol:** RPS6KA6

**Gene Alias:** RSK4