

Datasheet

RPS6KA6 monoclonal antibody (M04), clone 8E8

Catalog Number: H00027330-M04

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant RPS6KA6.

Clone Name: 8E8

Immunogen: RPS6KA6 (NP_055311.1, 636 a.a. ~ 745 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

RIGNGKFSLSGGNWDNISDGAKDLLSHMLHMDPHQR
YTAEQILKHSWITHRDQLPNDQPKRNDVSHVVKGAMV
ATYSALTHKTFQPVLEPVAASSLAQRRSMKKRTSTGL

Host: Mouse

Reactivity: Human, Rat

Applications: ELISA, IHC-P, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 27330

Gene Symbol: RPS6KA6

Gene Alias: RSK4