

Datasheet

MED31 monoclonal antibody (M01A), clone 2C8

Catalog Number: H00051003-M01A

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant MED31.

Clone Name: 2C8

Immunogen: MED31 (AAH12539, 1 a.a. ~ 131 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MAAAVAMETDDAGNRLRFQLELEFVQCLANPNYLNFL
AQRGYFKDKAFVNYLKYLLYWKDPEYAKYLKYPQCLH
MLELLQYEHFRKELVNAQCAKFIDEQQILHWQHYSRK
RMRLQQALAEQQQNNTSGK

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51003

Gene Symbol: MED31

Gene Alias: 3110004H13Rik, CGI-125, FLJ27436, FLJ36714, Soh1