

## Datasheet

### SH3GLB1 monoclonal antibody (M01), clone 1B3-A5

**Gene Symbol:** SH3GLB1

**Gene Alias:** Bif-1, CGI-61, KIAA0491, dJ612B15.2

**Catalog Number:** H00051100-M01

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a full length recombinant SH3GLB1.

**Clone Name:** 1B3-A5

**Immunogen:** SH3GLB1 (AAH07455, 1 a.a. ~ 365 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

MNIMDFNVKKLAADAGTFLSRAVQFTEEKLGQAEKTE  
LDAHLENLLSKAECTKIWTEKIMKQTEVLLQPNNARIE  
EFVYEKLDKAPSRINPELLGQYMIDAGTEFGPGTAY  
GNALIKCGETQKRIGTADRELIQTSALNFLTPLRNFIEG  
DYKTIAKERKLLQNKRLDLDAKTRLKAKAAETRNS  
EQELRITQSEFDRQAEITRLLLEGISSTHAHHLRCLNDF  
VEAQMTYYAQCYYMLDLQKQLGSFSPNYLSNNNQ  
SVTPVPSVLPNAIGSSAMASTSGLVITSPSNLSDLKEC  
SGSRKARVLYDYDAANSTELSLLADEVITVFSVVGMD  
DWMGERGNQKGKVPITYLELLN

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG1 kappa

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 51100