

## Datasheet

### FLJ20643 monoclonal antibody (M06), clone 1E8

**Catalog Number:** H00055011-M06

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a full length recombinant FLJ20643.

**Clone Name:** 1E8

**Immunogen:** FLJ20643 (AAH01108, 1 a.a. ~ 290 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

MANPKLLGLELSEAEAGADSARFEELLLQASKELQQA  
QTTRPESTQIQPQPGFCIKTNSSEGKVFINICHSPSIPP  
PADVTEEEELQMLEEDQAGFRIPMSLGEPHAELDAKG  
QGCTAYDVAVNSDFYRRMQNSDFLRELVITAREGLE  
DKYNLQLNPEWRMMKNRPFMGSSISQQNIRSEQRPRI  
QELGDLYTPAPGRAESGPEKPHLNLWLEAPDLLLAEID  
LPKLDGALGLSLEIGENRLVMGGPQQLYHLDAYIPLQI  
NSHESKAAFHRKRKQLMVAMPLLLVPS

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG1 Kappa

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 55011

**Gene Symbol:** PIH1D1

**Gene Alias:** FLJ20643, NOP17