

Datasheet

SH2D4A monoclonal antibody (M01), clone 3G8

Catalog Number: H00063898-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant SH2D4A.

Clone Name: 3G8

Immunogen: SH2D4A (AAH14525, 239 a.a. ~ 338 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

KKANELLLSTGMPGSFLIRVSERIKGYALSYLESDGCK
HFLIDASADAYSFLGVDQLQHATLADLVEYHKEEPITSL
GKELLYPCGQQDQLPDYLELFE

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 63898

Gene Symbol: SH2D4A

Gene Alias: FLJ20967, SH2A