

Datasheet

RASGEF1C monoclonal antibody (M01), clone 3H8

Catalog Number: H00255426-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant RASGEF1C.

Clone Name: 3H8

Immunogen: RASGEF1C (NP_778232, 168 a.a. ~ 260 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

QGPEGLVGADKPISYRTKPPASIHRELLGVCSDPYTLA
QLLTHVELERLRHIGPEEFVQAFVNDPLASTKPCFSD
KTSNLEAYVKWFRNLCY

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 255426

Gene Symbol: RASGEF1C

Gene Alias: FLJ35841