

Datasheet

POU5F1 monoclonal antibody, clone 7E7

Gene Alias: MGC22487, OCT3, OCT4, OTF3, OTF4

Catalog Number: MAB10378

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against synthetic peptide of POU5F1.

Clone Name: 7E7

Immunogen: A synthetic peptide corresponding to human POU5F1.

Host: Mouse

Theoretical MW (kDa): 45

Reactivity: Human

Applications: ELISA, Flow Cyt, IF, WB-Ce
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Isotype: IgG1

Recommend Usage: ELISA (1:10000)
Western Blot (1:500-1:2000)
Immunofluorescence (1:200-1:1000)
Flow cytometry (1:200-1:400)
The optimal working dilution should be determined by the end user.

Storage Buffer: In ascites (0.03% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5460

Gene Symbol: POU5F1