

Datasheet

Anti-Cetuximab monoclonal antibody, clone 6D4

Catalog Number: MAB11129

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against F(ab)₂ fragment of Cetuximab.

Clone Name: 6D4-9

Immunogen: F(ab)₂ fragment of Cetuximab.

Host: Mouse

Applications: Blocking, ELISA
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody blocks Cetuximab binding to A-431 (EGFR-expressing) cells in flow cytometry assay.

Form: Liquid

Purification: Protein G purification

Isotype: IgG1, kappa

Recommend Usage: ELISA (2 ug/mL for coating)
The optimal working dilution should be determined by the end user.

Storage Buffer: In TBS, pH 8.0 (0.09% sodium azide)

Storage Instruction: Store at 4°C for at least 12 months. For long term storage store at -80°C.
Aliquot to avoid repeated freezing and thawing.