

## Datasheet

### Anti-Rituximab monoclonal antibody, clone 10C5

**Catalog Number:** MAB11131

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against F(ab)2 fragment of Rituximab.

**Clone Name:** 10C5-3

**Immunogen:** F(ab)2 fragment of Rituximab.

**Host:** Mouse

**Applications:** Blocking, ELISA  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody is bivalent. It blocks the Rituximab-CD20 interaction in ELISA.

**Form:** Liquid

**Purification:** Protein G purification

**Isotype:** IgG1, kappa

**Recommend Usage:** ELISA (2 ug/mL for coating)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In TBS, pH 8.0 (0.09% sodium azide)

**Storage Instruction:** Store at 4°C for at least 12 months. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.