

Datasheet

Mdc1 monoclonal antibody, clone P2B11 (PerCP)

mKIAA0170

Catalog Number: MAB11475**Regulatory Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody raised against recombinant partial Mdc1.**Clone Name:** P2B11**Immunogen:** Recombinant GST fusion protein corresponding to C-terminus of mouse Mdc1.**Host:** Mouse**Reactivity:** Bovine, Chimpanzee, Human, Mouse**Applications:** IF, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Liquid**Conjugation:** PerCP**Purification:** Protein G purification**Isotype:** IgG1**Recommend Usage:** Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)**Storage Instruction:** Store at 4°C.
Aliquot to avoid repeated freezing and thawing.**Entrez GeneID:** 240087**Gene Symbol:** Mdc1**Gene Alias:** 6820401C03, AA413496, Nfbd1,