

Datasheet

Slc17a7 monoclonal antibody, clone S28-9 (PerCP)

Gene Symbol: Slc17a7

Gene Alias: BNPI, Vglut1

Catalog Number: MAB11527

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against recombinant of Slc17a7.

Clone Name: S28-9

Immunogen: Recombinant fusion protein corresponding to amino acids 493-560 at C-terminus of rat Slc17a7.

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Detects ~52 kDa. No cross-reactivity against Slc17a6.

Form: Liquid

Conjugation: PerCP

Purification: Protein G Purified

Isotype: IgG1

Recommend Usage: Immunohistochemistry

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol)

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 116638