

Datasheet

Slc17a6 monoclonal antibody, clone S29-29 (FITC)

Catalog Number: MAB11542

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against recombinant of Slc17a6.

Clone Name: S29-29

Immunogen: Recombinant fusion protein corresponding to amino acids 501-582 at C-terminus of rat Slc17a6.

Host: Mouse

Applications: IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Detects ~60 kDa.

Form: Liquid

Conjugation: FITC

Purification: Protein G Purified

Isotype: IgG1

Recommend Usage: Immunohistochemistry

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol)

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84487

Gene Symbol: Slc17a6

Gene Alias: Dnpi, Vglut2