

## Datasheet

### Slc17a6 monoclonal antibody, clone S29-29 (PerCP)

**Catalog Number:** MAB11545

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against recombinant of Slc17a6.

**Clone Name:** S29-29

**Immunogen:** Recombinant fusion protein corresponding to amino acids 501-582 at C-terminus of rat Slc17a6.

**Host:** Mouse

**Applications:** IHC, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Detects ~60 kDa.

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Protein G Purified

**Isotype:** IgG1

**Recommend Usage:** Immunohistochemistry

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol)

**Storage Instruction:** Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 84487

**Gene Symbol:** Slc17a6

**Gene Alias:** Dnpi, Vglut2