

Datasheet

Cacna1g monoclonal antibody, clone S178A-9 (PerCP)

Catalog Number: MAB11563

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against recombinant of Cacna1g.

Clone Name: S178A-9

Immunogen: Recombinant fusion protien corresponding to amino acids 2052-2172 at C-terminus of mouse Cacna1g.

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols