

Datasheet

Grin1 monoclonal antibody, clone S308-48 (PerCP)

Gene Symbol: Grin1

Gene Alias: NMDAR1, NR1

Catalog Number: MAB11725

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial synthetic protein of Grin1.

Clone Name: S308-48

Immunogen: Recombinant fusion protein corresponding to amino acids 42-361 at extracellular N-terminus of rat Grin1.

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein G purification

Isotype: IgG1

Recommend Usage: Immunohistochemistry (0.1-1 ug/mL)

Western Blot (1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS pH7.4 (50% glycerol and 0.09% sodium azide).

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 24408