

Datasheet

BrdU monoclonal antibody, clone SPM537

Catalog Number: MAB12073

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against Bromodeoxyuridine.

Clone Name: SPM537

Immunogen: Bromodeoxyuridine (BrdU) conjugated to KLH.

Host: Mouse

Theoretical MW (kDa): 307

Applications: IHC-P, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Protein A/G purification

Isotype: IgG1

Recommend Usage: Immunohistochemistry (0.5-1 ug/mL)

Western Blot (0.5-1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.05% BSA, 0.05% sodium azide).

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.