

Datasheet

Cd3e monoclonal antibody, clone 145-2C11 (PerCP/Cy5.5)

Catalog Number: MAB12205

Regulatory Status: For research use only (RUO)

Product Description: Armenian Hamster monoclonal antibody raised against Cd3e.

Clone Name: 145-2C11

Immunogen: Cd3e

Host: Hamster

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP/Cy5.5

Purification: Affinity Chromatography

Isotype: IgG

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

Entrez GeneID: 12501

Gene Symbol: Cd3e

Gene Alias: AI504783, CD3, CD3epsilon, T3e