

Datasheet

CD4 monoclonal antibody, clone RPA-T4 (FITC)

Catalog Number: MAB12208

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against CD4.

Clone Name: RPA-T4

Immunogen: CD4

Host: Mouse

Reactivity: Human

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Affinity Chromatography

Isotype: IgG1, kappa

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

Entrez GeneID: 920

Gene Symbol: CD4

Gene Alias: CD4mut