

## Datasheet

### **Cd4 monoclonal antibody, clone RM4-5 (PE/Cy7)**

**Catalog Number:** MAB12210

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against Cd4.

**Clone Name:** RM4-5

**Immunogen:** Cd4

**Host:** Rat

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PE/Cy7

**Purification:** Affinity Chromatography

**Isotype:** IgG2a, kappa

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

**Entrez GeneID:** 12504

**Gene Symbol:** Cd4

**Gene Alias:** L3T4, Ly-4