

Datasheet

Cd8a monoclonal antibody, clone 2.43 (PerCP/Cy5.5)

Catalog Number: MAB12214

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against Cd8a.

Clone Name: 2.43

Immunogen: Cd8a

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP/Cy5.5

Purification: Affinity Chromatography

Isotype: IgG2b

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

Entrez GeneID: 12525

Gene Symbol: Cd8a

Gene Alias: BB154331, Ly-2, Ly-35, Ly-B, Lyt-2