

Datasheet

Itgam monoclonal antibody, clone M1/70 (PE/Cy5)

Catalog Number: MAB12216

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against Itgam.

Clone Name: M1/70

Immunogen: Itgam

Host: Rat

Reactivity: Human, Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PE/Cy5

Purification: Affinity Chromatography

Isotype: IgG2b, kappa

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

Entrez GeneID: 16409

Gene Symbol: Itgam

Gene Alias: CD11b/CD18, CR3, CR3A, Cd11b, F730045J24Rik, Ly-40, MAC1, Mac-1, Mac-1a