

## Datasheet

### Itgam monoclonal antibody, clone M1/70 (PerCP/Cy5.5)

**Catalog Number:** MAB12217

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against Itgam.

**Clone Name:** M1/70

**Immunogen:** Itgam

**Host:** Rat

**Reactivity:** Human, Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP/Cy5.5

**Purification:** Affinity Chromatography

**Isotype:** IgG2b, kappa

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

**Entrez GeneID:** 16409

**Gene Symbol:** Itgam

**Gene Alias:** CD11b/CD18, CR3, CR3A, Cd11b, F730045J24Rik, Ly-40, MAC1, Mac-1, Mac-1a