

Datasheet

CD45.1 monoclonal antibody, clone A20 (FITC)

Gene Alias: B220, CD45R, Cd45, Ly-5, Lyt-4, T200, loc

Catalog Number: MAB12227

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against CD45.1.

Clone Name: A20

Immunogen: CD45.1

Host: Mouse

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: The clone A20 monoclonal antibody reacts with mouse CD45.1, it does not cross-react with CD45.2.

Form: Liquid

Conjugation: FITC

Purification: Affinity Chromatography

Isotype: IgG2a, kappa

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

Entrez GeneID: 19264

Gene Symbol: Ptprc