

## Datasheet

### **Ptprc monoclonal antibody, clone RA3-6B2 (APC)**

**Catalog Number:** MAB12228

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against Ptprc.

**Clone Name:** RA3-6B2

**Immunogen:** Abelson murine leukemia virus-induced pre-B tumor cells.

**Host:** Rat

**Reactivity:** Human, Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** APC

**Purification:** Affinity Chromatography

**Isotype:** IgG2a, kappa

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

**Entrez GeneID:** 19264

**Gene Symbol:** Ptprc

**Gene Alias:** B220, CD45R, Cd45, Ly-5, Lyt-4, T200, loc