

Datasheet

Ptprc monoclonal antibody, clone RA3-6B2 (PE)

Catalog Number: MAB12229

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against Ptprc.

Clone Name: RA3-6B2

Immunogen: Abelson murine leukemia virus-induced pre-B tumor cells.

Host: Rat

Reactivity: Human, Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PE

Purification: Affinity Chromatography

Isotype: IgG2a, kappa

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

Entrez GeneID: 19264

Gene Symbol: Ptprc

Gene Alias: B220, CD45R, Cd45, Ly-5, Lyt-4, T200, loc