

Datasheet

Ly6g monoclonal antibody, clone RB6-8C5 (PE/Cy7)

Catalog Number: MAB12236

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against Ly6g.

Clone Name: RB6-8C5

Immunogen: Native purified from normal murine bone marrow cells.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PE/Cy7

Purification: Affinity Chromatography

Isotype: IgG2b, kappa

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

Entrez GeneID: 546644

Gene Symbol: Ly6g

Gene Alias: Gr-1, Gr1, Ly-6G