

## Datasheet

### Ly76 monoclonal antibody, clone TER-119 (PE)

**Catalog Number:** MAB12237

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against Ly76.

**Clone Name:** TER-119

**Immunogen:** Ly76

**Host:** Rat

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PE

**Purification:** Affinity Chromatography

**Isotype:** IgG2b, kappa

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH7.2 (0.09% Sodium azide, may contain carrier protein/stabilizer).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.

**Entrez GeneID:** 104231

**Gene Symbol:** Ly76

**Gene Alias:** TER-119, Ter119