

Datasheet

NANOG monoclonal antibody, clone 60CT77.1.1

Catalog Number: MAB12279

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against human NANOG.

Clone Name: 60CT77.1.1

Immunogen: Recombinant protein corresponding to
human NANOG.

Host: Mouse

Reactivity: Human

Applications: IHC-P, WB-Ce
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Purification: Protein G purification

Isotype: IgG1, Kappa

Recommend Usage: Immunohistochemistry (1:10-1:50)
Western Blot (1:100-1:250)
The optimal working dilution should be determined by
the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store at 4°C. For long term
storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 79923

Gene Symbol: NANOG

Gene Alias: -