

Datasheet

CD4 monoclonal antibody, clone 4AT4 (PE-Cy5)

Catalog Number: MAB12375

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human CD4.

Clone Name: 4AT4

Immunogen: CD4

Host: Mouse

Reactivity: Human

Applications: Flow Cyt, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PE-Cy5

Purification: Protein A/G affinity chromatography

Isotype: IgG2b

Recommend Usage: Flow cytometry
Immunofluorescence
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide, 0.2% BSA)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 920

Gene Symbol: CD4

Gene Alias: CD4mut