

Datasheet

CD14 monoclonal antibody, clone 4AW4 (FITC)

Catalog Number: MAB12390

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human CD14.

Clone Name: 4AW4

Immunogen: CD14

Host: Mouse

Reactivity: Human

Applications: Flow Cyt, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein A/G affinity chromatography

Isotype: IgG1

Recommend Usage: Flow cytometry
Immunofluorescence
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide, 0.2% BSA)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 929

Gene Symbol: CD14

Gene Alias: -

Gene Summary: CD14 is a surface protein

preferentially expressed on monocytes/macrophages. It binds lipopolysaccharide binding protein and recently has been shown to bind apoptotic cells. Alternative splicing results in multiple transcript variants encoding the same isoform. [provided by RefSeq]