

Datasheet

NCAM1 monoclonal antibody, clone COC56 (FITC)

Catalog Number: MAB12409

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against human NCAM1.

Clone Name: COC56

Immunogen: NCAM1

Host: Mouse

Reactivity: Human

Applications: Flow Cyt, IF
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein A/G affinity chromatography

Isotype: IgG1

Recommend Usage: Flow cytometry
Immunofluorescence
The optimal working dilution should be determined by
the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide,
0.2% BSA)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 4684

Gene Symbol: NCAM1

Gene Alias: CD56, MSK39, NCAM