

## Datasheet

### NCAM1 monoclonal antibody, clone COC56 (PE)

**Catalog Number:** MAB12410

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against human NCAM1.

**Clone Name:** COC56

**Immunogen:** NCAM1

**Host:** Mouse

**Reactivity:** Human

**Applications:** Flow Cyt, IF  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PE

**Purification:** Protein A/G affinity chromatography

**Isotype:** IgG1

**Recommend Usage:** Flow cytometry  
Immunofluorescence  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (0.09% sodium azide, 0.2% BSA)

**Storage Instruction:** Store at 4°C. Do not freeze.

**Entrez GeneID:** 4684

**Gene Symbol:** NCAM1

**Gene Alias:** CD56, MSK39, NCAM