

Datasheet

FCGR1A monoclonal antibody, clone CT101 (FITC)

Catalog Number: MAB12414

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human FCGR1A.

Clone Name: CT101

Immunogen: FCGR1A

Host: Mouse

Reactivity: Human

Applications: Flow Cyt, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein A/G affinity chromatography

Isotype: IgG1

Recommend Usage: Flow cytometry
Immunofluorescence
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide, 0.2% BSA)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 2209

Gene Symbol: FCGR1A

Gene Alias: CD64, CD64A, FCRI, FLJ18345, IGFR1

Gene Summary: This gene encodes a protein that plays

an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1. [provided by RefSeq]