

Datasheet

Cd8a monoclonal antibody, clone 4AM43 (FITC)

Catalog Number: MAB12419

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against mouse Cd8a.

Clone Name: 4AM43

Immunogen: Cd8a

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein A/G affinity chromatography

Isotype: IgG2b

Recommend Usage: Flow cytometry
Immunofluorescence
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide, 0.2% BSA)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 12525

Gene Symbol: Cd8a

Gene Alias: BB154331, Ly-2, Ly-35, Ly-B, Lyt-2