

Datasheet

THY1 monoclonal antibody, clone 5E10

Gene Alias: CD90, FLJ33325

Catalog Number: MAB12871

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human THY1.

Clone Name: 5E10

Immunogen: Human erythroleukemia cell line (HEL).

Host: Mouse

Reactivity: Human

Applications: Flow Cyt, ICC, IHC-Fr, IHC-P, IP, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Protein A purification

Isotype: IgG1

Recommend Usage: Flow Cytometry
Immunocytochemistry
Immunohistochemistry (Frozen sections)
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (20 ug/mL)
Immunoprecipitation
Western Blot
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4.

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 7070

Gene Symbol: THY1