

## Datasheet

### c-Myc tag monoclonal antibody

**Catalog Number:** MAB1296

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against synthetic peptides of c-Myc tag.

**Immunogen:** A synthetic peptide (conjugated with KLH) corresponding to amino acids 408-432 of c-Myc tag.

**Host:** Mouse

**Applications:** WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Isotype:** IgG1

**Recommend Usage:** Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), 150 mM NaCl, pH 7.4 (50% glycerol, 0.02% sodium azide)

**Storage Instruction:** Store at 4°C. Do not freeze.