

Datasheet

Bromodeoxyuridine (BrdU) monoclonal antibody, clone MoBu-1

Catalog Number: MAB13234

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against Bromodeoxyuridine (BrdU).

Clone Name: MoBu-1

Immunogen: Bromodeoxyuridine (BrdU) conjugated with hemocyanine.

Host: Mouse

Applications: Flow Cyt, IF, IHC-P
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Protein A/G purification

Isotype: IgG1

Recommend Usage: Flow Cytometry (0.5-1 ug/million cells in 0.1 mL)

Immunofluorescence (0.5-1 ug/mL)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (0.5-1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS.

Storage Instruction: Store at -20 to -80°C.
Aliquot to avoid repeated freezing and thawing.