

## Datasheet

### **Bromodeoxyuridine (BrdU) monoclonal antibody, clone 85-2C8**

**Catalog Number:** MAB13236

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against Bromodeoxyuridine (BrdU).

**Clone Name:** 85-2C8

**Immunogen:** Bromodeoxyuridine (BrdU) conjugated with hemocyanine.

**Host:** Mouse

**Applications:** Flow Cyt, IF, IHC-P  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Protein A/G purification

**Isotype:** IgG1

**Recommend Usage:** Flow Cytometry (0.5-1 ug/million cells in 0.1 mL)

Immunofluorescence (0.5-1 ug/mL)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (0.5-1 ug/mL)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 10 mM PBS.

**Storage Instruction:** Store at -20 to -80°C.  
Aliquot to avoid repeated freezing and thawing.