

Datasheet

Thy1 monoclonal antibody, clone T1-A3

Catalog Number: MAB1356

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native Thy1.

Clone Name: T1-A3

Immunogen: Native from Con A-activated, B cell-depleted mouse LN cells.

Host: Mouse

Reactivity: Mouse

Applications: Flow Cyt, IP
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: React with mouse Thy-1.1 (allele-specific).

Form: Lyophilized

Isotype: IgG

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from PBS

Storage Instruction: Store at -20°C on dry atmosphere. After reconstitution with deionized water, store at -20°C or lower.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 21838

Gene Symbol: Thy1

Gene Alias: CD90, T25, Thy-1, Thy-1.2, Thy1.1, Thy1.2