

Datasheet

CD2 monoclonal antibody, clone TP1/31 (PerCP)

Gene Alias: SRBC, T11

Catalog Number: MAB13674

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human CD2.

Clone Name: TP1/31

Immunogen: Normal human peripheral blood lymphocytes.

Host: Mouse

Theoretical MW (kDa): 50

Reactivity: Human

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein A/G purification

Purity: > 90%

Isotype: IgG1

Recommend Usage: Flow Cytometry (20 μ L/10⁶ cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C.

Entrez GeneID: 914

Gene Symbol: CD2