

Datasheet

NCAM1 monoclonal antibody, clone B-A19 (FITC)

Gene Alias: CD56, MSK39, NCAM

Catalog Number: MAB13713

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human NCAM1.

Clone Name: B-A19

Immunogen: U937 cell line.

Host: Mouse

Theoretical MW (kDa): 175-185

Reactivity: Human

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein A/G purification

Purity: >90%

Isotype: IgG1

Recommend Usage: Flow Cytometry (20 μ L/ 10^6 cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C.

Entrez GeneID: 4684

Gene Symbol: NCAM1