

Datasheet

CD22 monoclonal antibody, clone HIB22 (FITC)

Gene Alias: FLJ22814, MGC130020, SIGLEC-2,
SIGLEC2

Catalog Number: MAB13846

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against human CD22.

Clone Name: HIB22

Immunogen: Tonsil cells.

Host: Mouse

Theoretical MW (kDa): 130-140

Reactivity: Human

Applications: Flow Cyt
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein A/G purification

Purity: >90%

Isotype: IgG1, kappa

Recommend Usage: Flow Cytometry (20 μ L/ 10^6 cells)
The optimal working dilution should be determined by
the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer,
0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid
prolonged exposure to light.

Entrez GeneID: 933

Gene Symbol: CD22