

Datasheet

CD33 monoclonal antibody, clone HIM3-4 (PerCP)

Gene Alias: FLJ00391, SIGLEC-3, SIGLEC3, p67

Catalog Number: MAB13878

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human CD33.

Clone Name: HIM3-4

Immunogen: KG1a cell line.

Host: Mouse

Theoretical MW (kDa): 67

Reactivity: Human

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Protein A/G purification

Purity: >90%

Isotype: IgG1, kappa

Recommend Usage: Flow Cytometry (20 μ L/ 10^6 cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.

Entrez GeneID: 945

Gene Symbol: CD33